Support arm – Operation and maintenance manual

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**en** The guarantee covers faults or defects in material or manufacture within a period of 5 years. Products are subject to minor technical modifications and design deviations. E. & O.E.

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# Content

Mounting instruction.....	3
Operational signage.....	3
Symbols used in this manual.....	4
General safety.....	5
Versions.....	6
Intended use.....	8
Intended user profile.....	8
Operating environment.....	8
Operating instructions.....	9
All support arms.....	9
Height adjustable support arms.....	9
Sideways adjustable support arms.....	10
Height and sideways adjustable support arms.....	10
Dismounting from wall track.....	11
Mounting on wall track.....	11
Dismounting from vertical track.....	12
Mounting on vertical track.....	12
Cleaning.....	13
To remove limescale.....	13
Maintenance.....	14
Trouble-shooting.....	14
Technical data.....	14
Dimensions, mm.....	15
Dimensions, inches.....	17
Product weights.....	19
Disposal and recycling.....	20

en

## Mounting instruction

A mounting instruction is included with the product

## Operational signage

An operational signage can be read from the operation instruction

# Symbols used in this manual

The following symbols are used in this manual as well as the associated documentation supplied with the product

WARNING indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury

CAUTION indicates a potentially hazardous situation which, if not avoided, may result in minor or moderate injury. It may also be used to alert against unsafe practices

IMPORTANT indicates important information about handling and use of the product

NOTICE is used to address practices not related to personal safety

Hand injury warning

# General safety

## **WARNING**

- Read and understand the manual completely before use Keep the manual close to the unit
- Comply with this manual to avoid accidents involving serious personal injury
- Never use or handle this product in other ways than are specified in this manual, as this can put your safety at risk and cause damage to the product.
- Make sure that everyone who is to install or use the product and its systems has the necessary information and access to this manual

## **CAUTION**

- Please comply with these instructions to avoid the product becoming damaged or destroyed.
- Do not perform repairs, disassembly, or assembly operations, extensions, re-adjustments or modifications to this product. These must be carried out by persons authorized by manufacturer only.
- Never allow children to operate this product unless under supervision The product is not intended for play
- Never use the product if it is defective or damaged. Do not use a defective product before it has been repaired

## **NOTICE**

- The information in this manual is based on correct fitting of the product in accordance with our Assembly instruction
- We cannot be held liable if the product is used in any way which differs from that described in the guide .
- We reserve the right to amend this manual without prior notice

# Versions

HA-L650

HB-L680


HC-L680

HD-L680

XHA-U622

XHB-U658

XHC-U658

XHD-U658

en

# Versions

FS-600  
FS-800  
FS-1000  
FS-1200

XFS-600  
XFS-800  
XFS-1000  
XFS-1200

FSL-900

# Intended use

The support arm is used as support for a person when he or she sits down or gets up from a toilet

## Intended user

The product can be used by all - adults, children, young and old. Children should be under the supervision of an adult when using this product.

## Operating environment

For indoor use in a wet environment, such as a bathroom or wash room, in private homes, nursing homes and hospitals. For institutional use, the expected life time is 3 years.

*Life cycle test and weight loading test performed according to ISO 17966:2016*

## Transport conditions

The product must be placed and stored in accordance with the packaging instructions. Store in a cool, ventilated and dry place.

### NOTICE

If the product is installed in an atmosphere heavy in salt or chloride, the warranty will be reduced from the normal 3 years to 1 year

Salt and chloride have a corrosive effect on the products components and will cause a decrease in the life expectancy of the product. The damaging effects can, however, be reduced if the product is showered with tap water on a daily basis.

### CAUTION


The surface of the product may become hot and / or cold if exposed to external sources of heat or cold (eg sunlight or frost)


# Operating instructions

## All support arms


Fold the support arm fully down. Fold it back up when not in use. You can hear a click when the support arm is correctly folded up

## en Height adjustable support arms


## Sideways adjustable support arms

Release the locking button ① and move the support arm sideways. Tighten the locking button again to secure the support arm in place

## Height and sideways adjustable support arms


en

Use the locking button ② and the locking button ③ to adjust the support arm horizontally and vertically

### Dismounting from wall track

To remove the support arm from the wall track, release the locking button ① and lift the support arm upwards


### Mounting on wall track


Hook the upper part of the support arm on the wall track

Click the lower part of the support arm firmly to the wall track Lock the locking button.



# Dismounting from vertical track

Remove the set screw ① and slide the support arm upwards and free of the vertical track

# Mounting on vertical track

Pinch the locking button and hold way which differs from that described in the guide arm down on the vertical track

Release the locking to secure the support arm in place .  
Make sure that the locking pin is engaged in the notch

Remount the set screw

# Cleaning

The product is to be cleaned with hot water or a mild soapy cleaning agent.

## CAUTION

Always read the material safety data sheet (MSDS) and supplied instructions of the cleaning, disinfecting or de-scaling agent and comply with instructions for safe use, dilution, rinsing, drying, handling and storage of the agent.

## IMPORTANT

- Do not use an autoclave or other type of steam cleaner.
- Avoid aggressive cleaning agents containing acids
- Never spray the product directly with a high pressure cleaner.
- Avoid use of abrasive, corrosive or chlorine-based cleaning or disinfecting products.

### To remove limescale

To remove limescale marks on the aluminium parts, use 3 parts 30% household vinegar to 7 parts water and afterwards wipe the parts with a cloth wrung out in clean water.

In countries with relatively high levels of lime in the water supply, it may be necessary to remove limescale deposits from the aluminium parts using a de-scaling agent such as acetic acid. Test this first on a area that is not visible when the product is in use.

# Maintenance

Remove dust from the wall mount of the sideways adjustable support arm once a month .Wipe with lubricant.

Check whether the connection part of the product is loose, whether there are cracks on the guide rail, whether the walltrack notch is severely worn, etc. Once these situations occur, you should stop using it and contact a professional to evaluate and repair it in time.

Expected service life of the support arm is 3 years

# Trouble-shooting



Symptom	Cause	Action
The product can't be moved sideways.	There might be dirt or dust on the horizontal walltrack.	Clean the walltrack and grease it with silicone.
The locking button doesn't lock	The lock is not properly adjusted.	See mounting instruction.



# Technical data



Max user weight	150 kg / 331 lbs
Materials	Anodized aluminium Stainless steel PA6 ABS



Product weight loading
The weight loading capacity of the product is directly related to the installed wall. Make sure that the wall is strong enough to bear a weight of at least 150 kg. Each mounting screw will bear the weight evenly.

Dimensions, mm

		
	HA-L650	XHA-U622
A	650	622
B	100	
C	210	
D	152	159

		
	HB-L680	XHB-U658
A	682	657
B	124	
C	424	
D	130	


		
	HC-L680	XHC-U658
A	682	657
B	104	
C	164	

		
	HD-L680	XHD-U658
A	721	697
B	124	
C	164	
D	424	


en

Dimensions, inches


	FS-600	FS-800
A	600	800
	FS-1000	FS-1200
A	1000	1200

	XFS-600	XFS-800
A	600-750	800-950
	XFS-1000	XFS-1200
A	1000-1150	1200-1350

en

	FSL-900
A	780
B	370
C	450
D	450

# Product weights

Model number	Product mass [kg]	Product mass (lbs)
HA-L650	3.50	7.71
HB-L680	5.30	11.68
HC-L680	4.00	8.81
HD-L680	6.00	13.22
XHA-U622	3.65	8.04
XHB-U658	5.90	13.00
XHC-U658	4.00	8.81
XHD-U658	6.50	14.33
FS-600	1.02	2.24
FS-800	1.32	2.91
FS-1000	1.60	3.52
FS-1200	1.86	4.10
XFS-600	1.22	2.68
XFS-800	1.55	3.41
XFS-1000	1.88	4.14
XFS-1200	2.21	4.87
FSL-900	1.50	3.30

en

# Disposal and recycling

The unit contains reusable materials. There are no known hazards or risks associated with disposal of the product. All components can be discarded after having been cleaned and disinfected

When disposing of the product, we recommend that it be disassembled and broken down into different waste groups for recycling or combustion

The main metal part of the product is a flame-resistant material, and the plastic part is a non-flammable material.

The following table provides information on the recycling and handling of the product.

Aluminium	Aluminium is theoretically 100% recyclable without any loss of its natural qualities . It requires only 5% of the energy used to produce primary aluminium. The quality of recycled aluminium is so high that it can be used again and again, even in the case of anodised products as the anodised layer is a natural part of the aluminium .
Steel	Steel can be recycled by being re-melted and included as a secondary material in the production of new steel .
TPE	Use plastic sorting and recycling systems where they are locally available, otherwise incinerating plastic has the added benefit of generating energy, which can be used for, for example, combined heat and power production.
ABS	As virgin ABS is somewhat expensive, recycling ABS is economically very attractive. .Recycled ABS can be blended with virgin material to produce products with lower cost while preserving the high quality

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